

# YOU TECARX

## YOUR FULLY CUSTOMIZABLE TECAR THERAPY

[www.youtecarx.com](http://www.youtecarx.com)



# CHOOSE THE IDEAL CONFIGURATION FOR YOUR NEEDS AND THOSE OF YOUR PATIENTS

## WHAT IS TECAR THERAPY? WHAT BENEFIT DOES IT OFFER?

CONTACT DIATHERMY – more commonly known as TECAR THERAPY – is a medical practice that applies radiofrequency to biological tissue with a frequency ranging from 300 kHz to 1 MHz without generating muscle contraction.

The position of the electrodes on the body determines the current path, thereby ensuring a deep, targeted, comfortable application that is also easy to control and manage.

The heat-induced vasodilation increases the exchange of substances promoting the drainage of the inflamed area and the influx of reparative cells.

Tecartherapy is based on the use of an electromedical device that uses high-frequency currents for the **CHEMICAL | THERMAL | MECHANICAL** biostimulation of tissues with different resistivity.

Physiotherapy, rehabilitation, and pain management are the fields that most use the beneficial effects of tecartherapy.

Tecar Therapy units are often the first piece of equipment to be bought.

THE NUMBER OF UNITS IS GROWING ALL OVER THE WORLD!

## WHAT FIELD DO YOU WORK IN?

- Sports Rehabilitation
- Dermal-Functional Rehabilitation
- Physio-Aesthetics
- Pelvic Floor Rehabilitation
- MuskuloSkeletal Rehabilitation:
  - Chronic & Acute
  - Post-Traumatic
  - Post-Surgical

## WHAT TECHNIQUES DO YOU COMBINE?

- Fascia Therapy
- Manual Therapy
- Therapeutic Exercise
- Osteopathy



# CONFIGURATION STEP BY STEP

## CREATE YOUR OWN TECAR THERAPY DEVICE

### BASE CONFIGURATION:


Power: 150 Watt  
Frequency: 500 KHz  
Output: Unipolar  
Emission: Continuous

### INCLUDES:

- 1 power cable
- 1 cable for unipolar applications
- 1 unipolar handpiece
- 1 steel neutral electrode
- 1 resistive electrode Ø 61 mm
- 1 capacitive electrode Ø 61 mm
- Cable with clamp for disposable neutral electrode
- Disposable gel adhesive electrode (5 pcs)



### CUSTOMIZATION OPTIONS



# THE PRICE? YOU CHOOSE!



BASE CONFIGURATION

CUSTOMIZATIONS

TOP MODEL

€

€€

€€€

€€€€

YOULECAR



PAGE 4



# CORE & SYSTEM

Choose the power, frequencies, and forms of emission



## CHOOSE THE POWER

150 WATT (included)

200 WATT

250 WATT



## CHOOSE THE FREQUENCIES

500 KHZ (included)

STANDARD PACK 3 FREQ. (500 KHz / 750 KHz / 1 MHz)

ON DEMAND PACK (500 KHz + a choice of three frequencies from 300 kHz to 1 MHz)



## CHOOSE THE EMISSION

CONTINUOUS (included)

SENSITIVE ATHERMY

HEA - HIGH ENERGY ATHERMY

ATHERMY PACK (SENSITIVE + HEA)



# AUX OUTPUT

The doorway to a whole world of new treatment possibilities

## BOOST THE POWER OF YOUTECARX!




The multifunction AUX output opens the door to the highest performance levels using special accessories and application systems.

**JOULE MODE** is an innovative operating method that includes a dedicated hardware and software system



# ACCESSORIES

Applicators designed to reach specific therapeutic targets





# UNIPOLAR

Complete your accessory kit for traditional unipolar mode operation  
(active electrode + neutral plate)






# AUTOMATIC APPS

Not operator-dependent



# NEUTRO DINAMICO®


Dynamic applications with two active handpieces and without neutral plate, using NeutroDinamico® accessories





# BIPOLAR

Work without neutral plate using  
Bipolar accessories




# INTIMITY

Pelvic Floor rehabilitation





# TRIGGER & FASCIA




# PHYSIO-AESTHETICS

Physio-aesthetics Accessories



# OCTOPUS CONTOUR

Extreme freedom for every kind of treatment





# PRACTICALITY



FOR MORE INFORMATION  
CONTACT OUR SPECIALISTS



By **Youmed** s.r.l - Via Enrico Pallini, 9 - Roma (Italy)

[www.youtecarx.com](http://www.youtecarx.com)

Youtecar X is a certified medical device

CE 1370